# CHIMIA REPORT/COMPANY NEWS

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Companies present themselves and their products

Please contact: Kretz AG, Postfach 105, CH-8706 Feldmeilen, info@kretzag.ch

## A Unique Filtration Device For The Removal Of Proteins And Lipids From Plasma Samples

### Abstract

Bioanalytical method development is often a trade-off between achieving a needed sensitivity and the desire for minimal sample analysis time. When

sensitivity is not paramount, methods based on dilution or protein precipitation are selected to reduce cost and analysis time. These methods, however, may not be as robust or clean as


more involved techniques such as solid phase extraction. We investigate the possibility of creating robust methods across samples and/or species by selective removal of matrix interferences including proteins, phospholipids, and lysophospholipids. This filtration based device, known as Captiva ND<sup>Lipids</sup>, gives improved cleanliness over protein precipitation, yet avoids additional method development or sample processing time. The method also normalizes disparate samples and may lead to more robust methods for the entire bioanalytical process.

### Determination of Clozapine and Norclozapine in Plasma

A mix of clozapine and norclozapine was spiked in plasma

at 32 ng/mL, filtered through the Captiva ND<sup>Lipids</sup> filtration plate, and run on a Pursuit<sup>TM</sup> XRs Ultra 2.8 C18 HPLC column.

The Captiva ND<sup>Lipids</sup> method yielded higher average absolute recoveries compared to samples processed by protein precipitation. Also, the fluctuation in recoveries over the range tested was much wider for protein-precipitated samples (e.g. norclozapine had recoveries ranging from 12.1 to 101.7), while filtered samples had a narrower range (41.6 to 76.7). This could be attributed to a possible higher plasma protein-binding capacity to both these analytes. The process of filtration through Captiva ND<sup>Lipids</sup> could possibly shift the equilibrium of bound plasma protein-drug complex

to the unbound state resulting in better recoveries and reproducibility.

### Universality of Method for Broad Ranging Analytes and Matrices

The unique membrane approach to lipid removal means that functional group tolerance to analytes within a plasma matrix is very good. Only when an analyte is of very high logP will performance reduce.

### Extended Column Lifetimes for Fast LC Applications

When combining Captiva ND<sup>Lipids</sup> with sub 2 µm or sub 3 µm columns, such as Pursuit™ XRs Ultra 2.8 and Pursuit UPS 2.4, back pressure issues become a thing of the past. A longevity study analyzing beta blockers, spiked in rat plasma after filtration with Captiva ND<sup>Lipids</sup>, shows virtually no changes in back pressure, retention time and peak shape with UV and MS detection. Even after more than 6000 injections the column characteristics remain the same. Using Captiva ND<sup>Lipids</sup> helps prolong column lifetime and therefore reduces downtime due to column changing and/or conditioning. The integrated combination of sub 2 µm or 3 µm columns with such a unique filtration device increases productivity and reduces cost, two crucial elements for high throughput laboratories. Captiva ND<sup>Lipids</sup> can act as a protective shield for analytical columns and detectors. Ion sources require less maintenance,

reducing disturbance of daily work flows.

### Simple Methodology and Ease of Use


Captiva ND<sup>Lipids</sup> ease of use lies in the use of a non-drip, lipid depleting filter. The non-drip feature allows for protein precipitation and protein removal in a single well. The lipid depletion provides cleanliness nearing that of traditional solid phase extraction (SPE). The removal of matrix interferences in the filtration mode also allows analyte concentration similar to SPE.

### Conclusions

- Captiva ND<sup>Lipids</sup> is a novel approach to a well adopted sample prep technique.
- The removal of phospholipids allows for reduced ion suppression and greater sensitivity.
- Relative recoveries are higher and more reproducible compared to standard precipitation methods.
- HPLC column lifetimes can be dramatically increased when adopting this technology.
- A broad range of analytes and plasma types can be used under very generic, facile conditions.

• Biopack Lab  
 Varian Distributor  
 Dr. Leroy Emmanuel  
 Av. Louis Ruchonnet 2  
 CH-1003 Lausanne  
 Tel. +41 21 213 03 15  
 Fax +41 21 544 13 64  
 eleroy@biopack.ch  
 www.biopack.ch

## Millennia Prime 532 nm CW DPSS Lasers

- 6, 10 and 15 W of CW output power at 532 nm
- Long life ProLite® diodes
- EternAlign™ for perfect alignment over life of laser
- Easy field serviceability
- Compact footprint to simplify integration
- Low noise of <0.04% rms
- Excellent pointing stability <2 µrad/°C

The Spectra-Physics Millennia® Prime™ Series is a family of rugged, high-power CW green DPSS lasers at 532 nm. With output power up to 15 W, Millennia Prime lasers are ideal for 24/7 applications where uptime is critical, including photovoltaic solar cell manufacturing, hard disk texturing, semiconductor inspection, material processing as well as scientific applications such as Ti:Sapphire, dye and solid-state laser pumping.

### Most Reliable

The new Millennia Prime lasers are designed specifically for 24/7 industrial use demanding high uptime. Each laser incorporates Spectra-Physics ProLite® diodes with industry leading reliability and lifetime. The laser optical resonator features EternAlign™ optical mount technology for perfect alignment over the lifetime of the laser. Key components (diode/fiber module) are field serviceable to shorten mean time to repair (MTTR). Millennia Prime lasers are tested against the toughest standards and back by a two year warranty, the best in the industry.

### Excellent Performance

The Millennia Prime lasers are available in 6, 10, 15 W versions. They feature exceptional low noise of <0.04% rms, TEM<sub>00</sub> output and excellent pointing stability at <2 µrad/°C. With a compact footprint, the lasers can be easily integrated into industrial machine tools or used where space is critical. The Millennia Prime series of lasers are CDRH and RoHS compliant.

With several thousands of units installed and in use today, the Millennia laser leads in product reliability, value and performance.

### OEM Applications

- Laser doping of solar cells
- Film subtitling
- Hard disk surface texturing
- Semiconductor wafer inspection
- Materials processing

### Scientific Applications

- Pumping CW and mode-locked Ti:Sapphire lasers
- Pumping solid state and dye lasers
- Spectroscopy
- Flow cytometry
- Medical diagnostics
- Photoluminescence
- Particle scattering
- Raman imaging

• Newport Spectra-Physics GmbH  
 Guerickeweg 7  
 D-64291 Darmstadt  
 Tel. +49-(0) 6151-708-0  
 Fax +49-(0) 6151-708-217  
 www.newport.com

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## NanoCalc-2000 Thin Film Measurement System

The NanoCalc-2000 Thin Film Measurement System is the ideal tool for fast, reliable and easy measurement of transparent and semi-transparent thin layers.

Depending on the layer and substrate material it can measure layers of just a few nanometers up to several hundreds of microns.

The large spectral range available from 250 nm (UV) to 1100 nm (NIR), as well as our special simulation algorithms allow the measurement not only of standard samples like oxides, nitrides and resists on ideal substrates like silicon wafer or glass but also the measurement of critical layers like DLC (diamond like carbon) on metal, layers on unpolished steel plates or the thickness of transparent foils.

Short measurement time and the very flexible positioning of the fiber head make this system the first choice for integration into existing machinery as well as for adaptation to almost any microscope and for in situ process control.

Measurement of multi-layers, 3D-thickness profiling with motorized 150 mm or 300 mm mapping stage, on-line measurement data acquisition and remote control via ActiveX are optionally available.

- GMP SA Laser & Photonics  
Dübendorfstrasse 11a  
CH-8117 Fällanden  
Tel. +41 44 825 34 00  
Fax +41 44 825 34 01  
www.gmp.ch

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- Metrohm Schweiz AG